

## Raccoonpox virus ACG19\_gp083 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX1011

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Raccoonpox virus ACG19\_gp083 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                  |
|----------------|------------------|
| Target Protein | ACG19_gp083      |
| Species        | Raccoonpox virus |
| Protein Class  | Virology         |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                              |
|-------------------|--------------------------------------------------------------|
| Target Protein    | ACG19_gp083                                                  |
| Full Name         | Entry-fusion complex component, myristylprotein              |
| Alternative Names | ACG19_gp083; Entry-fusion complex component, myristylprotein |
| Gene ID           | <a href="#">24528117</a>                                     |
| UniProt ID        | <a href="#">A0A0G3G4D9</a>                                   |