

Monkeypox virus D4L Expression Plasmid

Made to Order

Cat: MPYF-0722-KX1121

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus D4L expression plasmid. The product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Description

Target	D4L
Species	Monkeypox virus
Strain	Zaire-96-I-16
Vector	Expression vector or lentiviral vector
NCBI Reference	NC_003310.1
Bacterial Resistance	Based on vector
Mutations	Based on specific requirements
Form	Lyophilized
Storage	Ambient temperature

Target

Target	D4L
Full Name	D4L
Alternative Names	MPXVgp007; D4L
Gene ID	929049
UniProt ID	Q3I9K7