

Monkeypox virus E1R Expression Plasmid

Made to Order

Cat: MPYF-0722-KX1212

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus E1R expression plasmid. The product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Description

Target	E1R
Species	Monkeypox virus
Strain	Zaire-96-I-16
Vector	Expression vector or lentiviral vector
NCBI Reference	NC_003310.1
Bacterial Resistance	Based on vector
Mutations	Based on specific requirements
Form	Lyophilized
Storage	Ambient temperature

Target

Target	E1R
Full Name	E1R
Alternative Names	MPXVgp098; E1R; similar to Vaccinia virus strain Copenhagen D1R; mRNA capping enzyme, large subunit; RNA 5' triphosphatase and RNA guanylyl transferase activities
Gene ID	928883
UniProt ID	V9NSD6