

Monkeypox virus E6R Expression Plasmid

Made to Order

Cat: MPYF-0722-KX1217

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus E6R expression plasmid. The product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Description

Target	E6R
Species	Monkeypox virus
Strain	Zaire-96-I-16
Vector	Expression vector or lentiviral vector
NCBI Reference	NC_003310.1
Bacterial Resistance	Based on vector
Mutations	Based on specific requirements
Form	Lyophilized
Storage	Ambient temperature

Target

Target	E6R
Full Name	E6R
Alternative Names	MPXVgp103; E6R; similar to Vaccinia virus strain Copenhagen D6R; early transcript ion factor, VETF, small subunit
Gene ID	928933
UniProt ID	Q3I8R1