

Monkeypox virus A40L Expression Plasmid

Made to Order

Cat: MPYF-0722-KX1264

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus A40L expression plasmid. The product is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Description

Target	A40L
Species	Monkeypox virus
Strain	Zaire-96-I-16
Vector	Expression vector or lentiviral vector
NCBI Reference	NC_003310.1
Bacterial Resistance	Based on vector
Mutations	Based on specific requirements
Form	Lyophilized
Storage	Ambient temperature

Target

Target	A40L
Full Name	A40L
Alternative Names	MPXVgp150; A40L; similar to Vaccinia virus strain Copenhagen A38L; integral membrane glycoprotein; immunoglobulin-like; regulation of the influx of extracellular Ca ²⁺
Gene ID	928948
UniProt ID	Q3I8L0