

Monkeypox virus, B4R-F

Made to Order

Cat: MPYF-0722-KX1470

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus B4R Forward Primer

Product Description

Target	B4R
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	CDC Approved Sequences
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	B4R
Full Name	Hypothetical protein
Alternative Names	MPXVgp165; B4R; hypothetical protein; similar to vaccinia virus strain Copenhagen B2R and cowpox virus B2R; has ATP-binding domain and possible ATPase activity; the poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	928917
UniProt ID	V9NMD4