

## Monkeypox virus, D14L-R

Made to Order

Cat: MPYF-0722-KX1513

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus D14L Reverse Primer

### Product Description

|                            |                                      |
|----------------------------|--------------------------------------|
| Target                     | D14L                                 |
| Species                    | Monkeypox virus                      |
| Strain                     | Zaire-96-I-16                        |
| CDC/WHO Approved Sequences |                                      |
| Length (nt)                | 20 or based on specific requirements |
| Form                       | Lyophilized                          |
| Storage                    | Store at 2 - 8°C                     |

### Target

|                   |                                                                                                       |
|-------------------|-------------------------------------------------------------------------------------------------------|
| Target            | D14L                                                                                                  |
| Full Name         | D14L                                                                                                  |
| Alternative Names | MPXVgp017; D14L; similar to Vaccinia virus strain Copenhagen C3L; secreted complement binding protein |
| Gene ID           | <a href="#">928959</a>                                                                                |
| UniProt ID        | <a href="#">Q98VL5</a>                                                                                |