

Monkeypox virus, C7L-R

Made to Order

Cat: MPYF-0722-KX1529

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus C7L Reverse Primer

Product Description

Target	C7L
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	C7L
Full Name	Hypothetical protein
Alternative Names	MPXVgp033; C7L; hypothetical protein; similar to cowpox virus G1L and vaccinia virus strain Copenhagen F1L; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	928913
UniProt ID	Q3I8Y3