

Monkeypox virus, Q1L-R

Made to Order

Cat: MPYF-0722-KX1556

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus Q1L Reverse Primer

Product Description

Target	Q1L
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	Q1L
Full Name	Hypothetical protein
Alternative Names	MPXVgp060; Q1L; hypothetical protein; similar to vaccinia virus strain Copenhagen O1L; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a biphasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	929033
UniProt ID	V9NP95