

Monkeypox virus, G2L-R

Made to Order

Cat: MPYF-0722-KX1567

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus G2L Reverse Primer

Product Description

Target	G2L
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	G2L
Full Name	Hypothetical protein
Alternative Names	MPXVgp071; G2L; hypothetical protein; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins; the Entomopoxvirinae are a sub-family of the Poxviridae that infect insects
Gene ID	929062
UniProt ID	Q5IXV4