

Monkeypox virus, G10R-R

Made to Order

Cat: MPYF-0722-KX1575

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus G10R Reverse Primer

Product Description

Target	G10R
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	G10R
Full Name	Poxvirus myristoylprotein
Alternative Names	MPXVgp079; G10R; poxvirus myristoylprotein; virion protein covalently modified by myristic acid during translation; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	928969
UniProt ID	V9NPW6