

Monkeypox virus, L4R-R

Made to Order

Cat: MPYF-0722-KX1584

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus L4R Reverse Primer

Product Description

Target	L4R
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	L4R
Full Name	L4R
Alternative Names	MPXVgp088; L4R; similar to Vaccinia virus strain Copenhagen J4R; RNA polymerase, 22 kDa subunit
Gene ID	929064
UniProt ID	Q5IXT7