
Monkeypox virus, E3R-R

Made to Order

Cat: MPYF-0722-KX1596

This product is for research use only and is not intended for diagnostic use.

Monkeypox virus E3R Reverse Primer

Product Description

Target	E3R
Species	Monkeypox virus
Strain	Zaire-96-I-16
CDC/WHO Approved Sequences	CDC Approved Sequences
Length (nt)	20 or based on specific requirements
Form	Lyophilized
Storage	Store at 2 - 8°C

Target

Target	E3R
Full Name	E3R
Alternative Names	MPXVgp100; E3R; similar to Vaccinia virus strain Copenhagen D3R; virion core protein
Gene ID	929052
UniProt ID	V9NU08