

Monkeypox virus L1R Recombinant Protein Expressed in *E. coli*

Cat: MPYF-0722-KX17

This product is for research use only and is not intended for diagnostic use.

This product is a Monkeypox virus L1R recombinant protein expressed in *E. coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	L1R
Species	Monkeypox virus
Protein Class	Virology
Molecular Weight	19 kDa

Product Description

Expression Systems	<i>E. coli</i>
Protein Tag	His tag at the N-terminus
Protein Form	Soluble
Form	LYOPH
Reconstitution	Please see Certificate of Analysis for specific instructions.
Buffer	PBS, pH7.4
Application	Functional Assay; Protein-protein Interaction; ELISA; Western Blotting; etc.
Purity	>90% as determined by SDS-PAGE
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	L1R
Full Name	L1R
Alternative Names	MPXVgp085; L1R; similar to Vaccinia virus strain Copenhagen J1R
Gene ID	928938

