

Monkeypox virus (Zaire-96-I-16) L5L Recombinant Protein

Made to Order

Cat: MPYF-0722-KX192

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus (Zaire-96-I-16) L5L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	L5L
Species	Monkeypox virus
Strain	Zaire-96-I-16
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	L5L
Full Name	L5L
Alternative Names	MPXVgp089; L5L; similar to Vaccinia virus strain Copenhagen J5L; essential for virus multiplication
Gene ID	928930
UniProt ID	Q3T724