

## Monkeypox virus (Zaire-96-I-16) H4L Recombinant Protein

Made to Order

Cat: MPYF-0722-KX197

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus (Zaire-96-I-16) H4L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                 |
|----------------|-----------------|
| Target Protein | H4L             |
| Species        | Monkeypox virus |
| Strain         | Zaire-96-I-16   |
| Protein Class  | Virology        |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                                                                                        |
|-------------------|------------------------------------------------------------------------------------------------------------------------|
| Target Protein    | H4L                                                                                                                    |
| Full Name         | H4L                                                                                                                    |
| Alternative Names | MPXVgp094; H4L; similar to Vaccinia virus strain Copenhagen H4L; virion core RNA polymerase-associated protein, RAP 94 |
| Gene ID           | <a href="#">928896</a>                                                                                                 |
| UniProt ID        | <a href="#">V9NQG7</a>                                                                                                 |