

Monkeypox virus (Zaire-96-I-16) E13L Recombinant Protein

Made to Order

Cat: MPYF-0722-KX213

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus (Zaire-96-I-16) E13L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	E13L
Species	Monkeypox virus
Strain	Zaire-96-I-16
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	E13L
Full Name	E13L
Alternative Names	MPXVgp110; E13L; similar to Vaccinia virus strain Copenhagen D13L; protein needed for the formation of immature IMV; surface membrane
Gene ID	929047
UniProt ID	Q3I8Q4