

Monkeypox virus (Zaire-96-I-16) B17R Recombinant Protein

Made to Order

Cat: MPYF-0722-KX281

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Monkeypox virus (Zaire-96-I-16) B17R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	B17R
Species	Monkeypox virus
Strain	Zaire-96-I-16
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	B17R
Full Name	B17R
Alternative Names	MPXVgp178; B17R; ankyrin-like
Gene ID	929045
UniProt ID	V9NY05