

Camelpox virus (M-96) CamMLVgp165 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX297

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp165 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp165
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp165
Full Name	CMLV165
Alternative Names	CamMLVgp165; CMLV165; superoxide dismutase-like protein; similar to vaccinia virus strain Copenhagen A45R, variola major virus strain Bangladesh A48R
Gene ID	932496
UniProt ID	Q77FQ3