

Camelpox virus (M-96) CamMLVgp134 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX335

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp134 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp134
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp134
Full Name	CMLV134
Alternative Names	CamMLVgp134; CMLV134; putative myristoylated membrane protein; similar to vaccinia virus strain Copenhagen A16L, variola major virus strain Bangladesh A17L
Gene ID	932656
UniProt ID	Q8V2N9