

Camelpox virus (M-96) CamMLVgp099 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX359

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp099 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp099
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp099
Full Name	CMLV099
Alternative Names	CamMLVgp099; CMLV099; putative IMV envelope protein p35; similar to vaccinia virus strain Copenhagen H3L, variola major virus strain Bangladesh I3L
Gene ID	932622
UniProt ID	Q8V2S3