

Camelpox virus (M-96) CamMLVgp192 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX362

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp192 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp192
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp192
Full Name	CMLV192
Alternative Names	CamMLVgp192; CMLV192; similar to vaccinia virus strain Copenhagen B15R/B14R, variola major virus strain Bangladesh B13R and cowpox virus B13R
Gene ID	932618
UniProt ID	Q8V2I7