

Camelpox virus (M-96) CamMLVgp189 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX387

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp189 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp189
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp189
Full Name	Kelch-like protein
Alternative Names	CamMLVgp189; Kelch-like protein; Kelch motifs are 50-amino acid motifs represented by one beta-sheet blade; may be involved in remodelling and stabilization of the host cytoskeleton and host immune evasion; BTB-Kelch proteins have Kelch repeats that can interact with actin filaments; the poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding n

on-structural proteins involved in genome replication and late genes encoding the viral structural proteins

Gene ID	932577
---------	------------------------

UniProt ID	Q8V2J0
------------	------------------------
