

Camelpox virus (M-96) CamMLVgp205 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX401

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp205 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp205
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp205
Full Name	CMLV205
Alternative Names	CamMLVgp205; CMLV205; putative serpin I-1 protein; similar to vaccinia virus strain Copenhagen C12L, variola major virus strain Bangladesh B21R and cowpox virus B20R
Gene ID	932553
UniProt ID	Q8V2H4