

Camelpox virus (M-96) CamMLVgp066 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX478

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp066 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp066
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp066
Full Name	CMLV066
Alternative Names	CamMLVgp066; CMLV066; similar to vaccinia virus strain Copenhagen O1L, variola major virus strain Bangladesh Q1L
Gene ID	932568
UniProt ID	Q8V2V2