

Camelpox virus (M-96) CamMLVgp022 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX491

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Camelpox virus (M-96) CamMLVgp022 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CamMLVgp022
Species	Camelpox virus
Strain	M-96
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CamMLVgp022
Full Name	Hypothetical protein
Alternative Names	CamMLVgp022; hypothetical protein; similar to vaccinia virus strain copenhagen C4 L and cowpox virus C16L; The poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins
Gene ID	932520

