

Cowpox virus (Brighton Red) CPXV074 CDS Recombinant Protein

Made to Order

Cat: MPYF-0722-KX642

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Cowpox virus (Brighton Red) CPXV074 CDS recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CPXV074 CDS
Species	Cowpox virus
Strain	Brighton Red
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CPXV074 CDS
Full Name	CPXV074 protein
Alternative Names	CPXV074; CPXV074 CDS; CPXV074 protein
Gene ID	1485950
UniProt ID	Q8QMZ4