

Cowpox virus (Brighton Red) CPXV035 CDS Recombinant Protein

Made to Order

Cat: MPYF-0722-KX678

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Cowpox virus (Brighton Red) CPXV035 CDS recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CPXV035 CDS
Species	Cowpox virus
Strain	Brighton Red
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CPXV035 CDS
Full Name	CPXV035 protein
Alternative Names	CPXV035; CPXV035 CDS; CPXV035 protein
Gene ID	1485910
UniProt ID	Q8QN27