
Cowpox virus (Brighton Red) CPXV214 CDS Recombinant Protein

Made to Order

Cat: MPYF-0722-KX705

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Cowpox virus (Brighton Red) CPXV214 CDS recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	CPXV214 CDS
Species	Cowpox virus
Strain	Brighton Red
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	CPXV214 CDS
Full Name	CPXV214 protein
Alternative Names	CPXV214; CPXV214 CDS; CPXV214 protein
Gene ID	1486093
UniProt ID	Q7THV0