

Vaccinia virus (Western Reserve) VACWR156 Recombinant Protein Expressed in Yeast

Cat: MPYF-0722-KX73

This product is for research use only and is not intended for diagnostic use.

This product is a Vaccinia virus (Western Reserve) VACWR156 recombinant protein expressed in Yeast. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	VACWR156
Species	Vaccinia virus
Strain	Western Reserve
Protein Class	Virology
Molecular Weight	16.2 kDa

Product Description

Expression Systems	Yeast
Protein Tag	His tag at the N-terminus
Protein Length	ECD
Protein Form	Soluble
Form	LYOPH
Reconstitution	Please see Certificate of Analysis for specific instructions.
Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Purity	>85% as determined by SDS-PAGE
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	VACWR156
Full Name	EEV membrane phosphoglycoprotein
Alternative Names	VACWR156; A33R; EEV membrane phosphoglycoprotein

Gene ID [3707686](#)

UniProt ID [P68617](#)
