

Ectromelia virus (Moscow) EVM033 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX812

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Ectromelia virus (Moscow) EVM033 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	EVM033
Species	Ectromelia virus
Strain	Moscow
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	EVM033
Full Name	Ser/Thr kinase
Alternative Names	ECTVgp033; EVM033; Ser/Thr kinase
Gene ID	951565
UniProt ID	Q8JLH1