

## Ectromelia virus (Moscow) EVM141 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX832

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Ectromelia virus (Moscow) EVM141 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                  |
|----------------|------------------|
| Target Protein | EVM141           |
| Species        | Ectromelia virus |
| Strain         | Moscow           |
| Protein Class  | Virology         |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                        |
|-------------------|------------------------|
| Target Protein    | EVM141                 |
| Full Name         | EVM141                 |
| Alternative Names | ECTVgp142; EVM141      |
| Gene ID           | <a href="#">951539</a> |
| UniProt ID        | <a href="#">Q8JL78</a> |