

Ectromelia virus (Moscow) EVM005 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX858

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Ectromelia virus (Moscow) EVM005 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	EVM005
Species	Ectromelia virus
Strain	Moscow
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	EVM005
Full Name	Ankyrin repeat domain-containing protein
Alternative Names	ECTVgp005; EVM005; ankyrin repeat domain-containing protein
Gene ID	951511
UniProt ID	Q8JLJ3