

Ectromelia virus (Moscow) EVM037 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX866

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Ectromelia virus (Moscow) EVM037 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	EVM037
Species	Ectromelia virus
Strain	Moscow
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	EVM037
Full Name	EVM037
Alternative Names	ECTVgp037; EVM037
Gene ID	951499
UniProt ID	Q8JLG7