

Raccoonpox virus ACG19_gp001 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX899

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Raccoonpox virus ACG19_gp001 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	ACG19_gp001
Species	Raccoonpox virus
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	ACG19_gp001
Full Name	Chemokine binding protein
Alternative Names	ACG19_gp001; Chemokine binding protein
Gene ID	24528243
UniProt ID	A0A0G3FXF9