

Raccoonpox virus ACG19_gp151 Recombinant Protein

Made to Order

Cat: MPYF-0722-KX946

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Raccoonpox virus ACG19_gp151 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	ACG19_gp151
Species	Raccoonpox virus
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	ACG19_gp151
Full Name	C-type lectin-like IEV/EEV glycoprotein
Alternative Names	ACG19_gp151; C-type lectin-like IEV/EEV glycoprotein
Gene ID	24528185
UniProt ID	A0A0G3G075