

Anti-Monkeypox Virus G6R Monoclonal Antibody (MP-W30) (Mouse IgG)

Cat: MPYF-0822-YF65

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Monkeypox Virus G6R. The mAb with highest affinity will be selected. The selected antibody recognizes Monkeypox Virus G6R. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Warnings For Research Use Only. Not for Clinical Use or Diagnostic Use.

Product Description

Target	G6R
Species Reactivity	Monkeypox Virus
Strain	Strain Zaire-96-I-16
Specificity	This antibody recognizes monkeypox virus (MPXV) G6R.
Antibody Isotype	Mouse IgG
Clone	MP-W30
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	G6R
Gene Name	G6R
Introduction	Monkeypox virus (MPV) is a double-stranded DNA, zoonotic virus and a species of the genus Orthopoxvirus in the family Poxviridae. It causes a disease that has similar characteristics to smallpox but with a milder rash and a lower death rate. The virus causes disease in both humans and animals.
Alternative Names	Monkeypox virus; p28; D5R; E3 ubiquitin-protein ligase p28-like

Official Symbol	DNA-directed RNA polymerase 7 kDa subunit
-----------------	---

UniProt ID	Q77HN4
------------	------------------------