

---

## Anti-Camelpox virus CamMLVgp174 Monoclonal Antibody (MP-K151) (Mouse IgG)

**Cat: MPYF-0922-KX1016**

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Camelpox virus CamMLVgp174. The mAb with highest affinity will be selected. The selected antibody recognizes Camelpox virus CamMLVgp174. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

### Product Description

Target	CamMLVgp174
Species Reactivity	Camelpox virus
Strain	M-96
Cross Reactivity	CamMLVgp174
Specificity	This antibody recognizes Camelpox virus CamMLVgp174.
Antibody Isotype	Mouse IgG
Clone	MP-K151
Clonality	Monoclonal Antibody
Purification	≥95% (SDS-PAGE)
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

### Target

Target	CamMLVgp174
Gene Name	CamMLVgp174
Introduction	Camelpox is an economically important, contagious, often sporadic skin disease of camelids. The Camelpox virus (CMLV), which causes the disease, is closely linked to the Variola virus (VARV), which causes smallpox. Camelpox is confined to camel-rearing belts, primarily in developing nations, and has a negative economic impact due to significant illness, death, weight loss, and milk yield reductions.
Alternative Names	CamMLVgp174; Kelch-like protein; intracellular viral protein; may be involved in

remodelling and stabilization of the host cytoskeleton; similar to Vaccinia virus Copenhagen A55R; the poxviridae are enveloped unsegmented dsDNA viruses; unlike many dsDNA viruses that replicate in the host nucleus poxviruses encode their own replication machinery and therefore replicate in the cytoplasm; viral genes are expressed in a bi-phasic manner with early genes encoding non-structural proteins involved in genome replication and late genes encoding the viral structural proteins

---

Official Symbol	Truncated kelch-like protein
Gene ID	<a href="#">932685</a>
UniProt ID	<a href="#">Q8V2K4</a>

---