

Skunkpox virus (WA) BIZ96_gp088 Recombinant Protein

Made to Order

Cat: MPYF-0922-KX123

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Skunkpox virus (WA) BIZ96_gp088 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	BIZ96_gp088
Species	Skunkpox virus
Strain	WA
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	BIZ96_gp088
Full Name	crescent membrane and immature virion formation protein
Alternative Names	BIZ96_gp088; crescent membrane and immature virion formation protein; SKPV-WA-088
Gene ID	29057666
UniProt ID	A0A1C9KBM7