

## Skunkpox virus (WA) BIZ96\_gp066 Recombinant Protein

Made to Order

Cat: MPYF-0922-KX212

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Skunkpox virus (WA) BIZ96\_gp066 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

Target Protein	BIZ96_gp066
Species	Skunkpox virus
Strain	WA
Protein Class	Virology

### Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

### Target

Target Protein	BIZ96_gp066
Full Name	hypothetical protein
Alternative Names	BIZ96_gp066; hypothetical protein; SKPV-WA-066
Gene ID	<a href="#">29057644</a>
UniProt ID	<a href="#">A0A1C9KBL8</a>