

Taterapox Virus (Dahomey 1968) TATV_DAH68_131 Recombinant Protein

Made to Order

Cat: MPYF-0922-KX232

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Taterapox Virus (Dahomey 1968) TATV_DAH68_131 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	TATV_DAH68_131
Species	Taterapox Virus
Strain	Dahomey 1968
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	TATV_DAH68_131
Full Name	p4a precursor of IMV core protein 4a
Alternative Names	TATV_DAH68_131; p4a precursor of IMV core protein 4a
Gene ID	4238766
UniProt ID	Q0NP80