

## Vaccinia Virus (Smallpox Vaccine) (Western Reserve) A53R Recombinant Protein Made to Order

Cat: MPYF-0922-KX448

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Vaccinia Virus (Smallpox Vaccine) (Western Reserve) A53R recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |                                   |
|----------------|-----------------------------------|
| Target Protein | A53R                              |
| Species        | Vaccinia Virus (Smallpox Vaccine) |
| Strain         | Western Reserve                   |
| Protein Class  | Virology                          |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                    |
|-------------------|----------------------------------------------------|
| Target Protein    | A53R                                               |
| Full Name         | secreted TNF-receptor-like protein                 |
| Alternative Names | A53R; secreted TNF-receptor-like protein; VACWR179 |
| Gene ID           | <a href="#">3707710</a>                            |
| UniProt ID        | <a href="#">P24756</a>                             |