

Vaccinia Virus (Smallpox Vaccine) (Western Reserve) H3L Recombinant Protein Made to Order

Cat: MPYF-0922-KX533

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Vaccinia Virus (Smallpox Vaccine) (Western Reserve) H3L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	H3L
Species	Vaccinia Virus (Smallpox Vaccine)
Strain	Western Reserve
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	H3L
Full Name	temporal expression: late
Alternative Names	H3L; temporal expression: late; VACWR101
Gene ID	3707557
UniProt ID	A0A0M4RY97