

Vaccinia Virus (Smallpox Vaccine) (Western Reserve) G4L Recombinant Protein Made to Order

Cat: MPYF-0922-KX621

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Vaccinia Virus (Smallpox Vaccine) (Western Reserve) G4L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	G4L
Species	Vaccinia Virus (Smallpox Vaccine)
Strain	Western Reserve
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	G4L
Full Name	temporal expression: late
Alternative Names	G4L; temporal expression: late; VACWR081
Gene ID	3707537
UniProt ID	P68460