

Variola (Smallpox) Virus (India-1967, Ind3) K3L Recombinant Protein

Made to Order

Cat: MPYF-0922-KX695

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Variola (Smallpox) Virus (India-1967, Ind3) K3L recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	K3L
Species	Variola (Smallpox) Virus
Strain	India-1967, Ind3
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	K3L
Full Name	DNA-binding phosphoprotein
Alternative Names	K3L; DNA-binding phosphoprotein; VARVgp057
Gene ID	1486458
UniProt ID	P0DOP3