

Volepox virus (CA) BI079_gp040 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX152

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Volepox virus (CA) BI079_gp040 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	BI079_gp040
Species	Volepox virus
Strain	CA
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	BI079_gp040
Full Name	ribonucleotide reductase small subunit
Alternative Names	BI079_gp040; ribonucleotide reductase small subunit; VPXV-CA-040
Gene ID	29063987
UniProt ID	A0A1C9KC71