

Anti-Variola (Smallpox) Virus A29L Monoclonal Antibody (MP-K1567) (Mouse IgG)

Cat: MPYF-1022-KX1633

This product is for research use only and is not intended for diagnostic use.

The antibody customized production platform is used to screen the mouse monoclonal antibodies with at least 5 synthetic peptides of Variola (Smallpox) Virus A29L. The mAb with highest affinity will be selected. The selected antibody recognizes Variola (Smallpox) Virus A29L. The isotype is Mouse IgG. It can be used in applications: WB (Other applications need to be tested.).

Product Description

Target	A29L
Species Reactivity	Variola (Smallpox) Virus
Strain	India-1967, Ind3
Cross Reactivity	A29L
Specificity	This antibody recognizes Variola (Smallpox) Virus A29L.
Antibody Isotype	Mouse IgG
Clone	MP-K1567
Clonality	Monoclonal Antibody
Purity	≥95% (SDS-PAGE)
Purification	Purified by Protein A/G chromatography.
Applications	WB (Other applications need to be tested.)
Buffer	PBS, pH 7.4
Storage	Store at 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year). Avoid repeated freeze-thaw cycles.

Target

Target	A29L
Gene Name	A29L
Introduction	Variola virus is a sizable, brick-shaped virus that is between 302 and 350 nm and 244 and 270 nm in size. It has a single, linear, double-stranded DNA genome that is 186 kbp in size and has a hairpin loop at each end. Variola major and variola minor are the two traditional subtypes of smallpox. Although monkeys and other animals h

ave contracted the virus in an experimental environment, variola only infects humans in nature.

Alternative Names	A29L; Uncharacterized protein; VARVgp133
Official Symbol	Uncharacterized Protein
Gene ID	1486505
UniProt ID	Q89489