

Volepox virus (CA) BI079_gp166 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX35

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Volepox virus (CA) BI079_gp166 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	BI079_gp166
Species	Volepox virus
Strain	CA
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	BI079_gp166
Full Name	inactive Cu-Zn superoxide dismutase-like virion protein
Alternative Names	BI079_gp166; inactive Cu-Zn superoxide dismutase-like virion protein; VPXV-CA-166
Gene ID	29064113
UniProt ID	A0A1C9KCK5