

## Volepox virus (CA) BI079\_gp150 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX47

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Volepox virus (CA) BI079\_gp150 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

|                |               |
|----------------|---------------|
| Target Protein | BI079_gp150   |
| Species        | Volepox virus |
| Strain         | CA            |
| Protein Class  | Virology      |

### Product Description

|                    |                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Expression Systems | Based on specific requirements                                                                                                 |
| Protein Tag        | Based on specific requirements                                                                                                 |
| Protein Length     | Full length                                                                                                                    |
| Protein Form       | Soluble                                                                                                                        |
| Storage            | Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles. |

### Target

|                   |                                                  |
|-------------------|--------------------------------------------------|
| Target Protein    | BI079_gp150                                      |
| Full Name         | NT transdomain protein                           |
| Alternative Names | BI079_gp150; NT transdomain protein; VPXV-CA-150 |
| Gene ID           | <a href="#">29064097</a>                         |
| UniProt ID        | <a href="#">A0A1C9KCI0</a>                       |