

Pseudocowpox Virus (VR634) PCPV_gp129 Recombinant Protein

Made to Order

Cat: MPYF-1022-KX471

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Pseudocowpox Virus (VR634) PCPV_gp129 recombinant protein. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Target Protein	PCPV_gp129
Species	Pseudocowpox Virus
Strain	VR634
Protein Class	Virology

Product Description

Expression Systems	Based on specific requirements
Protein Tag	Based on specific requirements
Protein Length	Full length
Protein Form	Soluble
Storage	Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein	PCPV_gp129
Full Name	Ankyrin/F-box protein
Alternative Names	PCPV_gp129; Ankyrin/F-box protein
Gene ID	8797330
UniProt ID	D3IZF1